

## RATINGS/PERFORMANCE

<b>Type Number</b>		MLW-2422C	MLW-2425C
<b>LCD Diagonal Size</b>		24 inch 518 mm (H) x 324 mm (V)	
<b>Pixels</b>		1920 (H) x 1200 (V) 16:10	
<b>Brightness</b>		400cd/m <sup>2</sup>	
<b>Contrast Ratio</b>		1000:1	
<b>Viewing Angle (L/R)</b>		178°	
<b>Input Video Signal</b>	<b>HD-SDI</b>	1ch (BNC incl. looping out) *Option: M-MSDI2 board	1ch (BNC incl. looping out) *2 <sup>nd</sup> input available as option: M-MSDI2
	<b>DVI</b>	1 ch (DVI-D)	1 ch (DVI-D) *2 <sup>nd</sup> input available as option: M-RGB3
	<b>Composite, Y/C</b>	1ch (BNC, Y/C incl. looping out), NTSC/PAL	
	<b>Analogue Component</b>	1ch (BNC incl. looping out), RBG/Y, Pb, Pr	1ch (BNC incl. looping out), RBG/Y, Pb, Pr *2 <sup>nd</sup> input available as option: M-RGB3
	<b>PC</b>	1ch (D-Sub 15pin), Analogue R, G, B, (VGA to WUXGA)	1ch (D-Sub 15pin), Analogue R, G, B, (VGA to WUXGA) *2 <sup>nd</sup> input available as option: M-RGB3
<b>Dual Display</b>		N/A	Picture in Picture, Side by Side
<b>Video Process Engine</b>		10 bit	
<b>Colour Space</b>		102 % (of NTSC)	
<b>Colour Gradinet</b>		10 bit (1 Billion colours)	
<b>Power Requirement</b>		AC 100V - 120V, 200 - 240V ±10% 50/60Hz	
<b>Weight / Dimensions</b>		Approx. 7.2 kg; W 583 x H 442 x D 109 mm	Approx. 7.4 kg; W 583 x H 442 x D 109 mm