

## RATINGS / PERFORMANCE

<b>Type Number</b>		MLW-2626C	MLW-2626C-3D
<b>LCD Diagonal Size</b>		26" 576 mm (H) x 324 mm (V)	
<b>Pixels</b>		1920 (H) x 1080 (V) FHD	
<b>Brightness</b>		450 cd/m <sup>2</sup>	350 cd/m <sup>2</sup>
<b>Contrast Ratio</b>		1400:1	
<b>Viewing Angle (L/R)</b>		178°	
<b>3D effect</b>		N/A	Passive 3D (Side by Side, Line by Line, Top and Bottom, Simultaneous) with Xpol® polarization filter
<b>Input Video Signal</b>	<b>3D</b>	N/A	2ch (BNC incl. looping out)
	<b>HD-SDI</b>	2ch (Mulch SDI)	1ch (BNC incl. looping out)
	<b>DVI</b>	2ch (DVI-D)	
	<b>Composite, Y/C</b>	1ch (BNC, Y/C incl. looping out), NTSC/PAL	
	<b>Analogue Component</b>	1ch (BNC incl. looping out)	
	<b>PC</b>	1ch (D-Sub 15 pin), Analogue R, G, B, (VGA to WUXGA)	
<b>Multi Input</b>		Picture in Picture, Side by Side, Swap display is available	
<b>Video Process Engine</b>		10 bit	
<b>Colour Space</b>		72 % (of NTSC)	68 % (of NTSC)
<b>Colour Gradinet</b>		10 bit x 3 (1.06 Billion colours)	
<b>Power Requirement</b>		DC 24V/3.3A	
<b>Weight / Dimensions</b>		Approx. 8 kg / W 640 x H 418 x D 82 mm	