

RATINGS / PERFORMANCE MLW-2750UHD	
LCD size	27-inch
Panel type	IPS
Colour Gamut	sRGB 115%
Colour Depth	10bit colour
Brightness	800 cd/m ²
Aspect Ratio	16:9
Resolution	3840 x 1080
View angle	Horizontal : 178°, Vertical : 178°
Contrast ratio	1000:1
Response time	14ms
4K HDMI (2.0)	1ch
4K DisplayPort (1.2)	1ch, with loop-through
HD/PC 3G/HD-SDI	1ch, with loop-through
HD/PC DVI	1ch, with loop-through
USB	RS-232C
Others	1ch, with loop-through
Cable Cover (Rear)	Available
Power Input	100-240V, 50/60Hz
Power Consumption	120W
PIP	Available
PBP	Available
Safety	FDA, MDD, Medical UL, cUL, CE, KC
Water, Dust proof	(Front) IP35, (Rear) IP32
Weight	7.7kg
Dimensions	654.4 (W) x 410.9 (H) x 58 (D) mm
Monitor Fan	Fanless
Standard Accessory	Instruction manual, Power cable (2pcs), AC/DC adaptor, DisplayPort cable, HDMI cable