

RATINGS/ PERFORMANCE

Type Number		MLW-3225C	MLW-3225C-3D
LCD Diagonal Size		32" 698 mm (H) x 393 mm (V)	
Pixels		1920 (H) x 1080 (V) FHD	
Brightness		500 cd/m ²	200 cd/m ²
Contrast Ratio		1300:1	
Viewing Angle (L/R)		178°	
3D effect		N/A	Passive 3D (Line by Line) with Xpol® polarization filter
Input Video Signal	3D	N/A	2ch (BNC incl. looping out)
	HD-SDI	1ch (BNC incl. looping out) *2nd input available as option: M-MSDI2	1ch (BNC incl. looping out)
	DVI	1ch (DVI-D) *2nd input available as option: M-RGB3	
	Composite, Y/C	1ch (BNC, Y/C incl. looping out), NTSC/PAL	
	Analogue Component	1ch (BNC incl. looping out), RBG/Y, Pb, Pr *2nd input available as option: M-RGB3	
	PC	1ch (D-Sub 15 pin), Analogue R, G, B, (VGA to WUXGA) *2nd input available as option: M-RGB3	
Multi Input		Picture in Picture, Side by Side	
Video Process Engine		10 bit	
Colour Space		72% (of NTSC)	
Colour Gradinet		10 bit (1 Billion colours)	
Power Requirement		AC 100V - 120V, 200 - 240V ±10% 50/60Hz	
Weight / Dimensions		Approx. 13 kg / W780 x H493 x D117.5 mm	